




Armin Baack

Supply for Laboratory and Medicin


The competent partner of your laboratory

Agagel Horizontal Gel Electrophoresis Apparatus



Delivery : -chamber
 -casting stand
 -2 combs
 -1 UV-transparent gel-tray
 combs and trays at buyer's option
 ready to use, leak-free

Agagel Horizontal Gel Electrophoresis Apparatus				
description	order number	gel size	buffer volume	price
Micro-horizontal	BAA26900240	7x4cm	180ml	290,00EUR/each
Mini-horizontal	BAA26900025	10,5x8,2cm	400ml	340,00EUR/each
Midi-horizontal	BAA26900012	14,5x12cm	550ml	360,00EUR/each
Midi-large-horizontal	BAA26900001	17x15,1cm	850ml	390,00EUR/each
Midi large Bufferflow	BAA26900090	17x15,1cm	850ml	430,00EUR/each
Maxi-horizontal	BAA26900200	25x20,4cm	1150ml	495,00EUR/each
Maxi-wide-horizontal	BAA26900260	14,5x23,9cm	1050ml	495,00EUR/each
Maxi-large-horizontal	BAA26900060	25x23,9cm	1350ml	580,00EUR/each
Maxi-xl-horizontal	BAA26900280	40x23,9cm	5400ml	680,00EUR/each



- System is ready to use without clamps
- Very good cost/performance ratio
- No leakage at all, easy Handling
- 4mm retractable plugs, gold coated
- Specially protected cables, savety lid
- Platinum electrodes

Accessories (chombs,...) you can find: www.baacklab.de

Speedy 1 Agagel Horizontal Mini Gel Apparatus

- for fast analysis of DNA and RNA probes
- preparing the gel directly in apparatus
- low consumption of agarose and buffer
- good running of gels
- best results by 80V, fast by 130V (1%)

BAA26900046 with 2 combs 320,00EUR/each





Armin Baack

Supply for Laboratory and Medicin

The competent partner of your laboratory

Polyacrylamide Vertical Gel Electrophoresis Apparatus



H 10 without cooling base



H10 with cooling base




H10 casting stand

H10 mini

- recommended for short and fast runs
- buffer volume 400-800ml
- weight incl. glass plates and casting stand
- 4mm retractable plugs, gold coated
- platinum electrodes
- safety lid and leakage free box
- rubber feeds

order- no	Description	EUR/each
BAA2691000	H10 Mini vertical, incl. 2pairs of glassplates, casting stand, 2 combs, 4 spacer on glass, buffer volume 400-800ml	590,-



H15 Midi vertikal

- recommended for medium runs
- buffer volume : 600 - 1200 ml
- weight incl. glass plates and casting stand
- 4 mm retractable plugs, gold


coated

- specially protected cables
- platinum electrodes
- safety lid and leakage free box
- rubber feeds

BAA26910050	H15 Midi vertical, incl. 2 pairs of glassplates, casting stand,	a.A,-
	2 combs, 4 spacer, buffer volume 600 - 1200 ml	

BAA26910100	H16 Midi wide, incl. 2 pairs of glassplates, casting stand,	a.A
	2 combs, 4 spacer, buffer volume 2000 - 4000 ml	

BAA26910150	H 20 Maxi vertical, incl. 2 pairs of glassplates, casting stand,	980,-
	2 combs, 4 spacer, buffer volume 2000 - 4000 ml	



H20 Maxi vertikal

- recommended for medium runs with many samples
- buffer volume : 2000 - 4000 ml
- weight incl. glass plates and casting stand
- 4 mm retractable plugs, gold

coated

- specially protected cables
- platinum electrodes
- safety lid and leakage free box
- rubber feeds

BAA26910200	Vario vertical, incl. 1 pair of glassplates, 1 comb	a.A.
	3 spacer, buffer volume 1000 - 1500 ml, Dm: 27x27x60cm (LxWxH)	



More information:

www.baacklab.de

prices 2021 +VAT+transport